

## Auto-Titrate Comparison Chart

	Auto-Titrate CPAP Machines		Auto-Titrate BiLevel Machines	
	S8 AutoSet II	REMstar Auto M Series	VPAP™ Auto 25	BiPAP Auto
<b>General</b>				
Manufacturer	ResMed	Respironics	ResMed	Respironics
Warranty	2 year	2 year	2 year	2 year
Noise Level	<24±2 dBA @ 10 cm	Not available	<26 dB	Not available
Weight	3.09 lbs	2.2 lbs	2.9 lbs	2.2 lbs
Dimensions	5.7"L x 6.5"W x 4.6"H	7.5"L x 5.0"W x 3.125"H	4.4"L x 6.5"W x 5.7"H	7.5"L x 5.0"W x 3.1"H
<b>Electrical</b>				
110/120 Volt	YES	YES	YES	YES
220/250 Volt	YES	YES	YES	YES
DC Cable	YES	YES	YES	YES
Battery Operation	w/converter	YES	YES	YES
Integrated Battery	NO	NO	NO	NO
<b>Functional Features</b>				
Heated Humidifier	H4i	M Series	H4i	M Series
Cold Humidifier	H4i	M Series	H4i	M Series
Pressure Range	4 to 20 cm	4 to 20 cm	3 to 25 cm	4 to 25 cm
Exhale Pressure Relief	YES - EPR	YES - C-Flex/A-Flex	YES - EPR	YES - BiFlex
Altitude Adjustment	Automatic	Automatic	Automatic	Automatic
Ramp	YES	YES	YES	YES
Hr/Session Recording	YES	YES	YES	YES
Optional Software	NO	YES	NO	YES
Data Card	YES	YES	YES	YES
Leak Compensation	YES	YES	YES	YES
Mask Off Alert	YES	YES	YES	YES
Auto ON/OFF	YES	YES	YES	YES