

## BiLevel Machine Comparison Chart

	BiPAP Plus M Series	VPAP S
<b>General</b>		
Manufacturer	Respironics	ResMed
Warranty	2 year	2 year
Noise Level	Not available	<26 dB
Weight	2.2 lbs	2.9 lbs
Dimensions	1.5"L x 5.0"W x 3.1"H	4.4"L x 6.5"W x 5.7"H
<b>Electrical</b>		
110/120 Volt	YES	YES
220/250 Volt	YES	YES
DC Cable	YES	YES
Battery Operation	YES	YES
Integrated Battery	NO	NO
<b>Functional Features</b>		
Heated Humidifier	M Series	H4i
Cold Humidifier	M Series	H4i
Pressure Range	4 to 25 cm	3 to 25 cm
Exhale Pressure Relief	NO	YES
Altitude Adjustment	Automatic	Automatic
Ramp	YES	YES
Hr/Session Recording	YES	YES
Optional Software	NO	NO
Data Card	NO	YES
Leak Compensation	YES	YES
Mask Off Alert	YES	YES
Auto ON/OFF	YES	YES